# **Antithrombin III Activity**

Test ID: 216 CPT: 85300

# **Clinical Significance:**

Diagnosis of antithrombin deficiency, acquired or congenital

Monitoring treatment of antithrombin deficiency disorders, including infusion of antithrombin therapeutic concentrate

#### **Profile Components:**



# Labtech Diagnostics

# **Container:**

3.2% sodium citrate (light bluetop) tube

#### Transport Temperature:

Frozen

#### Specimen:

frozen platelet-poor plasma

#### Specimen Stability:

Room temperature: Unacceptable Refrigerated: Unacceptable

Frozen: 30 days

## Reject Criteria:

Hemolysis • Received room temperature • Received refrigerated

## **Days Performed:**

Mon-Fri; Report available: 3 days

#### **Collection Instructions:**

Centrifuge light blue-top tube for 15 minutes at approximately 1500 g within 60 minutes of collection. Using a plastic pipette, remove plasma, taking care to avoid the WBC/platelet buffy layer and place into a plastic vial. Centrifuge a second time and transfer platelet-poor plasma into a new plastic vial. Plasma must be free of platelets (<10,000/uL). Freeze immediately and ship Froozen.