

Androstenedione, LC/MS/MS

Test ID: 17182

CPT: 82157

Clinical Significance:

Androstenedione is useful when evaluating patients with androgen excess and managing patients with Congenital Adrenal Hyperplasia (CAH).

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST)

Transport Temperature:

Refrigerated

Specimen:

Serum

Specimen Stability:

Room temperature: 1 days
Refrigerated: 1 week
Frozen: 6 Months

Reject Criteria:

Serum Separator Tube (SST®) •
Glass tubes

Days Performed:

Sun-Fri

Collection Instructions:

Serum Separator Tube. 1 mL serum. Allow SST to clot in an upright position for at least 30 minutes, then centrifuge sample within 2 hours of collection. Clearly label tube as serum. Refrigerate. Sample should be collected between 6am and 10am.