

Alkaline Phosphatase, Bone Specific (BAP)

Test ID: 29498

CPT: 84075

Clinical Significance:

Quantitative measurement of skeletal alkaline phosphatase to be used as an aid in the management of patients with Paget disease and osteoporosis

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 7 days
Refrigerated: 7 days
Frozen: 60 days

Reject Criteria:

Gross hemolysis, Gross lipemia,
Non-serum samples, Grossly icteric

Days Performed:

Mon-Fri

Collection Instructions:

Allow serum samples to clot completely before centrifugation.