
Alpha-Globin Common Mutation Analysis

Test ID:	11175	СРТ:	81257
Clinical Significance:			
•	•	e among individuals of Asian heritage. Disease . Genetic counseling may be advised for some	



Container:	Transport Temperature:	Specimen:
EDTA (lavender-top) tube	Room temperature	Whole blood
Specimen Stability:	Reject Criteria:	Days Performed:
Room temperature: 7 days	Received frozen	Mon
Refrigerated: 7 days		

Collection Instructions:

Frozen: Unacceptable

Physician Attestation of Informed Consent

This germline genetic test requires physician attestation that patient consent has been received if ordering medical facility is located in AK, DE, FL, GA, IA, MA, MN, NV, NJ, NY, OR, SD or VT or test is performed in MA.