
Amylase, Random Urine with Creatinine

Test ID:	8464	CPT:	82150, 82570
Clinical Significance:			
	ul in the consideration of macroamylasem tration is elevated. With pseudocysts, the eatitis.		· ·



Container:	
Sterile screw-cap container	
·	

Transport Temperature:

Room temperature

Specimen:	
10 mL random urine	

Specimen Stability:
Room temperature: 7 days
Refrigerated: 7 days
Frozen: 28 days

Reject Criteria: Acidified urine Days Performed:

Mon-Fri

Collection Instructions:

Submit 10 mL of a well-mixed random collection. Collect without preservatives. Refrigerate during and after collection. Aliquot unpreserved specimen prior to any addition of acid.