

# Arsenic, Blood

Test ID: 269

CPT: 82175

## Clinical Significance:

This assay is used to monitor exposure to arsenic, wellness, and therapy during treatment of chronic myelocytic leukemia

## Profile Components:



# Labtech Diagnostics

### Container:

EDTA (royal blue-top) tube

### Transport Temperature:

Room temperature

### Specimen:

Whole blood

### Specimen Stability:

Room temperature: 10 days  
Refrigerated: 10 days  
Frozen: Unacceptable

### Reject Criteria:

Clotted specimen

### Days Performed:

Mon, Wed, Fri 8 a.m., Release 6 p.m.; Report available: 4 days

## Collection Instructions:

Avoid seafood consumption for 48 hours prior to sample collection.  
clean skin prior to venipuncture. Avoid worksite collection.

Carefully