

Autoimmune Hepatitis Diagnostic Panel

Test ID: 19873

CPT: 86038, 86255, 83516, 86376

Clinical Significance:

The detection of a number of autoantibodies supports the diagnosis of autoimmune hepatitis.

Profile Components:

ANA Screen (IFA) with Reflex to Titer, Mitochondrial Antibody with Reflex to Titer, Actin (Smooth Muscle) Antibody (IgG), Liver Kidney Microsome (LKM-1) Antibody (IgG)

If ANA Screen, IFA is positive, then ANA Titer and Pattern will be performed at an additional charge (CPT code(s): 86039).

If Mitochondrial Antibody Screen is positive, then Mitochondrial Antibody Titer will be performed at an additional charge (CPT code(s): 86256).



Container:

Serum Separator Tube (SST®) (x2)

Transport Temperature:

Refrigerated

Specimen:

Serum

Specimen Stability:

Room temperature: 4 hrs

Refrigerated: 7 days

Frozen: 30 days

Reject Criteria:

Received room temperature

Days Performed:

Mon-Fri

Collection Instructions: