

Allergen, Pinto Bean (f300) IgE

Test ID: 3083

CPT: 86003

Clinical Significance:

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 14 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Fri

Collection Instructions:

Separate serum from cells within 45 minutes of venipuncture.