

# Allergen, Turkey Meat (f284) IgE

Test ID: 2748

CPT: 86003

## Clinical Significance:

## Profile Components:



Labtech Diagnostics

### Container:

Serum Separator Tube (SST®)

### Transport Temperature:

Room temperature

### Specimen:

Serum

### Specimen Stability:

Room temperature: 14 days  
Refrigerated: 28 days  
Frozen: 1 year

### Reject Criteria:

### Days Performed:

Mon-Fri

## Collection Instructions:

Separate serum from cells within 45 minutes of venipuncture.