Antibody Screen, RBC, Reflex to Identification, Titer, and Antigen Typing

Test ID:	795	СРТ:	86850
Clinical Significance:			
	ct significant RBC antibodies.		

Profile Components:

If Antibody Screen is positive, Antibody Identification, Titer and Antigen Typing will be performed at an additional charge (CPT code(s): 86870 (x) # of panel(s) performed, 86886 (x) # of titer(s) performed and 86905 (x) # of antigen typing(s).Based on reactions presenting from Antibody Identification, then 36668 - Direct Antiglobulin Test (DAT) with Reflex to Anti-C3 and Anti-IgG may be performed at an additional charge (CPT code(s): 86880). **Direct Antiglobulin Test (DAT)** Is positive, then %36669 - Anti-C3 and Anti-IgG will be performed at an additional charge (CPT code(s): 86880 x2).

Container:

EDTA (lavender-top) tube

Transport Temperature:

Refrigerated

Specimen: whole Blood

Specimen Stability:

Room temperature: 24 hours Refrigerated: 7 days Frozen: Unacceptable

Reject Criteria:

Gross hemolysis • Grossly lipemic • Serum Separator Tube (SST[®]) • Received frozen • Cord blood • Grossly icteric **Days Performed:**

Mon-Fri

Collection Instructions: