
## Antiproteinase 3 (PR3) Abs

Test ID:	9918	CPT:	86021
Clinical Significance:			
Evaluating patients suspe	cted of having Wegener granuloma	atosis (WG)	
<b>U U</b>		conjunction with MPO / Myeloperoxidase Ar / be obtained as VASC / Antineutrophil Cytopl	
May be useful to follow tr	eatment response or to monitor d	lisease activity in patients with myeloperoxida	ase antibodies



Transport Temperature:

Room temperature

Specimen:			
Serum			

Specimen Stability:	
Room temperature: 7 days	
Refrigerated: 14 days	
Frozen: 30 days	

Reject Criteria: Gross hemolysis

Days	Performed:

Mon-Fri

## **Collection Instructions:**

Separate serum from cells within 45 minutes of venipuncture.