## Aldosterone, Serum

Test ID:	17181	CPT:	82088
Clinical Significance:			
Investigation of primary alo	Josteronism (eg, adrenal adenor	ma/carcinoma and adrenal cortical hyperplasia) a	ind secondary
aldosteronism (renovascular disease, salt depletion, potassium loading, cardiac failure with ascites, pregnancy, Bartter		icy, Bartter	
syndrome)			

Profile Components:	
	Labtech Diagnostics

Container:
Red-top tube (no gel) or
EDTA (lavender-top)

Transport Temperature:

Room temperature

Specimen:	
Serum or Plasma	

Specimen Stability:
Room temperature: 4 days
Refrigerated: 7 days
Frozen: 28 days

Reject Criteria:

Days Performed:
Mon-Fri

## **Collection Instructions:**

Separate serum or plasma from cells within 45 minutes of venipuncture and lable it.