

# AA- Allergy, Dairy Panel

Test ID: 99992

CPT: 82785 , 86003 x 8


## Clinical Significance:

See individual tests

## Profile Components:

Total, IgE  
Egg White, f1, IgE  
Milk, f2, IgE  
Egg Yolk, f75, IgE  
Alpha-Lactalbumin, f76, IgE

Beta-lactoglob, f77, IgE  
Casein, f78, IgE  
Cheese Cheddar, f81, IgE  
Yogurt, f360, IgE



Labtech Diagnostics

## Container:

x2 Serum Separator Tube (SST®)

## Transport Temperature:

Room Temperature

## Specimen:

Serum

## Specimen Stability:

Room temperature: 28 days  
Refrigerated: 28 days  
Frozen: 1 year

## Reject Criteria:

## Days Performed:

Mon-Sat

## Collection Instructions: