

AA- Allergy, Grains & Beans Panel

Test ID: 99991

CPT: 82785 , 86003 x 17


Clinical Significance:

See individual tests

Profile Components:

Barly, F6
Gluten, F79
Grain
Green Bean, F315
Oat, F7
Rice, F9
Rye, F5

Sesame, F10
Soybean, F14
Wheat, F4



Labtech Diagnostics

Container:

x2 Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: