

AA- Allergy, Seafood Panel

Test ID: 9948


CPT: 82785 , 86003 x 14

Clinical Significance:

See individual tests

Profile Components:

Codfish, F3
Crab, F23
Shrimp, F24
Tuna, F40
Salmon, F41
Lobster, F80
Oyster, F290
Sea Food



Labtech Diagnostics

Container:

x2 Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: