

AA- Allergy, Trees Panel

Test ID: 99996

CPT: 82785 , 86003 x 18

Clinical Significance:

See individual tests

Profile Components:

Cottonwood, T14
Mountain Cedar, T6
Oak, T7
Pecan Tree, T22
Red Cedar, T219
Sycamore, T11
Walnut, T10
Willow, T12



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: