

ABO Group and Rh Type

Test ID: 7788

CPT: 86900, 86901

Clinical Significance:

ABO type and Rh are needed to identify candidates for Rh immune globulin and to assess the risk of hemolytic disease of the newborn.

Profile Components:

ABO Group and Rh Type



Labtech Diagnostics

Container:

EDTA (lavender-top) tube

Transport Temperature:

Room temperature

Specimen:

5 mL whole blood collected in an EDTA (lavender-top) tube

Specimen Stability:

Room temperature: 1 day
Refrigerated: 1 week
Frozen: Unacceptable

Reject Criteria:

-Hemolysis
-Samples collected in SST

Days Performed:

Mon-Fri

Collection Instructions:

5 mL EDTA whole blood. Room temperature.

When multiple tests are ordered, submit separate tube for this test.