

Aldosterone, 24-Hour Urine

Test ID: 229

CPT: 82088

Clinical Significance:

Evaluation of hypertension and suspected primary or secondary aldosteronism, usually in conjunction with urine creatinine, sodium and/or potassium levels, as well as a plasma renin level.

Profile Components:

Includes:
Creatinine



Labtech Diagnostics

Container:

Urine collection container

Transport Temperature:

Frozen

Specimen:

Urine

Specimen Stability:

Room Temperature: Unacceptable
Refrigerated: Unacceptable
Frozen: 1 month

Reject Criteria:

-Room temperature samples
-Refrigerated samples

Days Performed:

Sun, Tues, Thurs

Collection Instructions:

24 hour urine refrigerated during collection. Record total volume. Transfer 10 mL from a well-mixed 24 hour collection to a plastic transport tube. Freeze. Critical Frozen (CFZ).

When multiple tests are ordered, submit separate tube for this test.