

Allergen, House dust, Holl (h2) IgE

Test ID: 3539

CPT: 86003

Clinical Significance:

Latin name: Hollister-Stier Labs.

Mite, cockroach and animal allergens are major components of house dust. For more information, see Insects and Mites.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: