

Allergen, Kiwi Fruit (f84) IgE

Test ID: 2884

CPT: 86003

Clinical Significance:

Latin name: Actinidia deliciosa

Source material: Fresh fruit

Family: Actinidiaceae (Dilleniaceae)

Common names: Kiwi, Chinese gooseberry, Kiwifruit, Monkey peach, Sheep peach

Synonyms: A. latifolia var. deliciosa, A. chinensis deliciosa

There are 2 common species of Kiwi commercially available:

A. chinensis - Gold kiwi

A. deliciosa - Green kiwi (Chinese gooseberry)

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: