

Allergen, Lamb (f88) IgE


Test ID: 2888

CPT: 86003

Clinical Significance:

Code: f88
Latin name: Ovis spp.
Source material: Raw meat
Mutton is regarded as a weak allergen and preparations of lamb meat have been suggested as an alternative to hydrolyzed cow's milk formulas.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: