

# Allergen, Mackerel (f206) IgE

Test ID: 2713

CPT: 86003

## Clinical Significance:

Latin name: *Scomber scombrus*

Source material: Fresh fish

Family: Scombridae

In a study of fish- and crustacea-allergic adults, the reactivity to mackerel was the second highest.

Allergen Exposure

Mackerel is usually eaten cooked, marinated or smoked.

## Profile Components:



# Labtech Diagnostics

### Container:

Serum Separator Tube (SST®)

### Transport Temperature:

Room Temperature

### Specimen:

Serum

### Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

### Reject Criteria:

### Days Performed:

Mon-Sat

## Collection Instructions: