

# Allergen, Malt (f90) IgE

Test ID: 2863

CPT: 86003

## Clinical Significance:

Latin name: Hordeum vulgare

Family: Poaceae (Gramineae)

Common names: Malt, Barley Malt

Food

A food, which may result in allergy symptoms in sensitised individuals.

## Profile Components:



# Labtech Diagnostics

### Container:

Serum Separator Tube (SST®)

### Transport Temperature:

Room Temperature

### Specimen:

Serum

### Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

### Reject Criteria:

### Days Performed:

Mon-Sat

## Collection Instructions: