

Allergen, Meadow fescue, (g4) IgE

Test ID: 3558

CPT: 86003

Clinical Significance:

Latin name: Festuca pratensis

Source material: Pollen

Family: Poaceae (Graminae)

Sub family: Pooideae

Tribe: Poeae

Common names: Meadow Fescue, Tall Fescue, Taller Fescue, Giant Fescue, English Bluegrass

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: