

Allergen, Melons (f87) IgE

Test ID: 2887

CPT: 86003

Clinical Significance:

Code: f87
Latin name: Cucumis melo spp.
Source material: Fruit
Family: Cucurbitaceae
Common names: Melon, Common melon, Muskmelon, Armenian cucumber

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: