

Allergen, Amoxicillin (c6) IgE


Test ID: 38476

CPT: 86003

Clinical Significance:

test is used for detection of IgE sensitization to the major allergenic determinant of amoxicillin, with potential cross-reactivity to other penicillins.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: