

Allergen, Mouse Epithelium (e71) IgE

Test ID: 2657

CPT: 86003

Clinical Significance:

Latin name: Mus spp.

Source material: Epithelium

Family: Muridae

Common names: Mouse, House mouse, Common house mouse

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: