

# Allergen, Oak, (t7) IgE

Test ID: 3094

CPT: 86003

## Clinical Significance:

Latin name: Quercus alba

Source material: Pollen

Family: Fagaceae

Common names: Oak, White oak, Forked-leaf white oak, Fork-leaf oak

## Profile Components:



# Labtech Diagnostics

### Container:

Serum Separator Tube (SST®)

### Transport Temperature:

Room Temperature

### Specimen:

Serum

### Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

### Reject Criteria:

### Days Performed:

Mon-Sat

## Collection Instructions: