

# Allergen, Oat, (f7) IgE

Test ID: 9338

CPT: 86003

## Clinical Significance:

Latin name: Avena sativa  
Source material: Untreated planting seeds  
Family: Graminae (Poaceae)  
Common names: Oats, Oat, Oatmeal, Oat groats

## Profile Components:



# Labtech Diagnostics

### Container:

Serum Separator Tube (SST®)

### Transport Temperature:

Room Temperature

### Specimen:

Serum

### Specimen Stability:

Room temperature: 28 days  
Refrigerated: 28 days  
Frozen: 1 year

### Reject Criteria:

### Days Performed:

Mon-Sat

## Collection Instructions: