
## Allergen, Apple (f49) IgE

Test ID:	2849	CPT:	86003
Clinical Significance:			
Distinguish between po	llen related apple allergy and fruit all	ergy due to LTP sensitization	
Improve the risk assess	ment using allergen components		
Improve management of	of apple allergic patients		



Conta	iner:			
Serum	Separat	tor Tube	e (SST®)	
	•		. ,	

Transport Temperature:

Room Temperature

Specimen:		
Serum		

Specimen Stability:		
Room temperature: 28 days		
Refrigerated: 28 days		
Frozen: 1 year		

Reject Criteria:

Days Performed:	
Mon-Sat	

**Collection Instructions:**