## Allergen, Rye Grass (g5) IgE

Test ID:	3100	СРТ:	86003
Clinical Significance:			
Latin name: Lolium perenne			
Source material: Pollen			
Family: Poaceae (Gramineae)			
Sub family: Pooideae			
Tribe: Poeae			
Common names: Rye grass, Ry	ye-grass, Perennial Rye Grass, P	erennial Rye-grass, Perennial Ryegrass, Ray-g	rass, Annual Ryegrass



Contai			-® )
Serum	Separato	r Tube (SST	°)

Transport Temperature:

Room Temperature

Specimen:		
Serum		

Specimen Stability:	
Room temperature: 28 days	
Refrigerated: 28 days	
Frozen: 1 year	

Reject Criteria:

Days Performed:	
Mon-Sat	

**Collection Instructions:**