

# Allergen, Asparagus (f261) IgE

Test ID: 2626

CPT: 86003

## Clinical Significance:

Asparagus is a widely grown vegetable. Together with Garlic, the Onion and the Leek, it belongs to the Alliaceae (previously known as Liliaceae) or Lily family, which contains ferns and vegetables and flowers such as tulips. The genus Asparagus is made up of some 300 species.

## Profile Components:



# Labtech Diagnostics

### Container:

Serum Separator Tube (SST®)

### Transport Temperature:

Room Temperature

### Specimen:

Serum

### Specimen Stability:

Room temperature: 28 days  
Refrigerated: 28 days  
Frozen: 1 year

### Reject Criteria:

### Days Performed:

Mon-Sat

## Collection Instructions: