

# Allergen, Plum (f255) IgE

Test ID: 3081

CPT: 86003

## Clinical Significance:

Latin name: Prunus domestica  
Source material: Fresh fruit  
Family: Rosaceae  
Common names: Plum, Gage, Prune  
Food  
A food, which may rarely result in allergy symptoms in sensitised individuals.

## Profile Components:



# Labtech Diagnostics

### Container:

Serum Separator Tube (SST®)

### Transport Temperature:

Room Temperature

### Specimen:

Serum

### Specimen Stability:

Room temperature: 28 days  
Refrigerated: 28 days  
Frozen: 1 year

### Reject Criteria:

### Days Performed:

Mon-Sat

## Collection Instructions: