

Allergen, Raspberry (f343) IgE

Test ID: 26281

CPT: 86003

Clinical Significance:

Latin name: *Rubus idaeus*

Source material: Fruit

Family: Rosaceae


Common names: Raspberry, Red raspberry, Wild raspberry, Common red raspberry, European red raspberry, American red raspberry

Synonyms: *R. buschii*, *R. vulgatus* var. *buschii*

Important included species:

R. occidentalis – Black raspberry/Thimbleberry

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: