

Allergen, Redtop(Bentgrass) (g9) IgE

Test ID: 3532

CPT: 86003

Clinical Significance:

Latin name: *Agrostis stolonifera*

Source material: Pollen

Family: Poaceae (Gramineae)

Sub family: Pooideae

Tribe: Agrostideae

Common names: Redtop, Bentgrass, Water Bent grass, Creeping Bent, Creeping Bentgrass, Carpet Bentgrass; also see below.

Synonyms: *A. lachnanthis*, *A. gigantea*, *A. alba*

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: