

Allergen, *Rhizopus nigricans* (m11) IgE

Test ID: 3734

CPT: 86003

Clinical Significance:

Latin name: *Rhizopus nigricans*

Source material: Spores and mycelium

Family: Mucoraceae


Mould

A mould, which may result in allergy symptoms in sensitised individuals.

Allergen Exposure

Fungi can be found throughout the world. They may live or be found in indoor as well as outdoor environments.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: