

Allergen, Rye, (f5) IgE

Test ID: 3556

CPT: 86003

Clinical Significance:

Latin name: Secale cereale

Source material: Untreated planting seeds

Family: Poaceae (Gramineae)

Common names: Rye, Rogge

There is a need to differentiate between Rye the foodstuff (Secale cereale), Cultivated Rye grass pollen (Secale cereale) g12, Rye grass pollen (Lolium perenne) g5, and Wild Rye grass pollen (Elymus tricoides) g70. This distantly related genus, Elymus, contains species known as Wild Ryes, which are used as cover and for forage. A Wheat hybrid is known as Giant Rye.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: