Allergen, Salmon (f41) IgE

Test ID:	2841	CPT:	86003
Clinical Significance:			
Latin name: Salmo salar			
Source material: Fish muscle			
Family: Salmonidae			
Confusion surrounds the com	non names of fish in this group	o. The genus Salmo salar is from the northern	Atlantic, the others
from the Pacific.			
Allergen Exposure			
Meat from salmon.			



Con	tair	ner:					
Ser	um :	Sepa	arato	r Tul	be (SS	ST®)	

Transport Temperature:

Room Temperature

Specimen:		
Serum		

Specimen Stability:
Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:
Mon-Sat

Collection Instructions: