

Allergen, Sardine/Pilchard (f61) IgE

Test ID: 2623

CPT: 86003

Clinical Significance:

Latin name: *Sardinops melanosticta*

Source material: Fish muscle

Family: Clupeidae

The name sardine or pilchard is given to herring-like fishes belonging to the genera *Sardina*, *Sardinella* and *Sardinops*. Canned sardines may, in fact, be herrings.

Allergen Exposure

The name sardine or pilchard is given to herring-like fishes belonging to the genera *Sardina*, *Sardinella* and *Sardinops*. Canned sardines may, in fact, be herrings.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: