

Allergen, Scallop (f338) IgE

Test ID: 9950

CPT: 86003

Clinical Significance:

Latin name: Pecten spp.

Source material: Muscle

Family: Pectinidae

Allergen Exposure

Scallops, or fan shells, are found in all seas, shallow as well as deep. Adult scallops are usually free-swimming. The well-developed muscle is the food served as scallop or, in France Coquille St. Jacques.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: