

Allergen, Shrimp, (f24) IgE

Test ID: 3016

CPT: 86003

Clinical Significance:

Latin name: *Pandalus borealis*, *Penaeus monodon*, *Metapenaeopsis barbata*, *Metapenaeus joyneri*
Source material: Boiled, frozen Atlantic shrimp and raw, frozen prawns from the Indo-West-Pacific
Family: Pandalidae (*Pandalus borealis*) and Penaeidae (*Penaeus monodon*, *Metapenaeopsis barbata*, *Metapenaeus joyneri*)
Shrimp is considered to be a highly allergenic food, causing severe reactions, such as anaphylaxis.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: