

Allergen, Sunflower Seed (k84) IgE

Test ID: 23864

CPT: 86003

Clinical Significance:

Latin name: Helianthus annuus

Source material: Whole seeds

Family: Asteraceae (Compositae)

Allergy to sunflower seeds is frequently associated with inhalant allergy to pollen from plants of the family Compositae.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: