

Allergen, Bahia grass, (g17) IgE

Test ID: 3048

CPT: 86003

Clinical Significance:

Common names: Bahia Grass, Bahiagrass, Paraguay Paspalum

A grass species producing pollen, which often induces hay fever, asthma and conjunctivitis in sensitised individuals

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: