

# Allergen, Thyme (f273) IgE

Test ID: 3076

CPT: 86003

## Clinical Significance:

Latin name: *Thymus vulgaris*

Family: Lamiaceae (Labiatae)

Common names: Thyme, Garden thyme, Common thyme


Herb

A herb, which may result in allergy symptoms in sensitised individuals.

Allergen Exposure

*Thymus* is a huge genus, containing 300 to 400 species, most of which are aromatic shrubs or perennials. All are native to Asia or Europe, most probably to Southern Europe. Thyme is now cultivated all over the world and has become naturalised in some areas, including the northeastern US.

## Profile Components:



# Labtech Diagnostics

### Container:

Serum Separator Tube (SST®)

### Transport Temperature:

Room Temperature

### Specimen:

Serum

### Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

### Reject Criteria:

### Days Performed:

Mon-Sat

## Collection Instructions: