

Allergen, Timothy Grass (g6) IgE

Test ID: 3072

CPT: 86003

Clinical Significance:

Latin name: Phleum pratense

Source material: Pollen

Common names: Timothy, Herd's Grass, Cat's Tail

Synonyms: P. nodosum, P. parnassicum

The two commonly recognised varieties are

P. Pratense var. pratense

P. Pratense var. nodosum

A grass species producing pollen, which often induces hayfever, asthma and conjunctivitis in sensitised individuals.

Profile Components:



Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: