

# Allergen, *Trichophyton rubrum* (m205) IgE

Test ID: 3406

CPT: 86003

## Clinical Significance:

Latin name: *Trichophyton* spp.

Source material: Spores and mycelium

Common names: *Trichophyton*

The following 3 allergens are included in this overview:

*Trichophyton rubrum* m205


*Trichophyton ment.* var *goetzii* m210

*Trichophyton ment.* var *interdigitale* m211

Mold

A mold, which may result in allergy symptoms in sensitised individuals.

## Profile Components:



# Labtech Diagnostics

### Container:

Serum Separator Tube (SST®)

### Transport Temperature:

Room Temperature

### Specimen:

Serum

### Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

### Reject Criteria:

### Days Performed:

Mon-Sat

## Collection Instructions: