
Allergen, Tuna (F40) IgE

Test ID:	2840	СРТ:	86003			
Clinical Significance:						
Latin name: Thunnus albacares	∟atin name: Thunnus albacares					
Source material: Fish muscle						
Family: Scombridae	Family: Scombridae					
Anaphylaxis after ingestion of tuna has been reported						
Allergen Exposure						
Meat from tuna.						
Tuna is sold fresh, canned or frozen.						
Japanese sashimi/sushi includes raw tu	าล.					



Container:			
Serum Sep	arator T	ube (SST®)
		·	

Transport Temperature:

Room Temperature

	Specimen Stability:			
Room temperature: 28 days				
	Refrigerated: 28 days			
	Frozen: 1 year			

Reject Criteria:

Days Performed:	
Mon-Sat	

Collection Instructions: