
Allergen, Whey (f236) IgE

Test ID:	3248	CPT:	86003
Clinical Significance:			
Latin name: Bos spp.			
Source material: Whey powe	der		
Whey is the by-product of cl	heese-making. Milk proteins cont	tain about 20% whey components.	
Allergen Exposure			
Processed cheese.			



Container:		
Serum Separator Tube (SST [®])		

Transport Temperature:

Room Temperature

Specimen:		
Serum		

Specimen Stability:		
Room temperature: 28 days		
Refrigerated: 28 days		
Frozen: 1 year		

Reject Criteria:

Days Performed:		
Mon-Sat		

Collection Instructions: