

# Allergen, White Ash (t15) IgE

Test ID: 2515

CPT: 86003

## Clinical Significance:

Latin name: Fraxinus americana

Source material: Pollen

Family: Oleaceae

There are 4 important genera in the Oleaceae family: Olive (Olea), Ash (Fraxinus), Lilac (Syringa), and Privet (Ligustrum).

## Profile Components:



# Labtech Diagnostics

### Container:

Serum Separator Tube (SST®)

### Transport Temperature:

Room Temperature

### Specimen:

Serum

### Specimen Stability:

Room temperature: 28 days

Refrigerated: 28 days

Frozen: 1 year

### Reject Criteria:

### Days Performed:

Mon-Sat

## Collection Instructions: