Allergen, White-faced hornet,(i2) IgE

| Test ID: | 3304 | CPT: | 86003 |
|---------------------------------|--------------------------------|---|--------------|
| | | | |
| Clinical Significance: | | | |
| Latin name: Dolichovespula ma | culata | | |
| Source material: Venom | | | |
| Family: Vespidae | | | |
| Order: Hymenoptera | | | |
| Some sting victims have been re | eported as suffering 2000 stin | ngs and surviving while others receive a single | fatal sting. |



| Container: | |
|--|--|
| Serum Separator Tube (SST [®]) | |
| | |
| | |
| | |
| | |
| | |

Transport Temperature:

Room Temperature

| Specimen: | |
|-----------|--|
| Serum | |
| | |
| | |
| | |
| | |

Reject Criteria:

| Days Performed: | |
|-----------------|--|
| Mon-Sat | |
| | |

Collection Instructions: